

**Please note:**

A location = Keizerzaal & Foyer

B location = Pické 1 room

 = BOOTH AVAILABLE (8m2)

 = BOOTH AVAILABLE (6m2)

 = BOOTH RESERVED


 = BOOTH ON OPTION

**Please note:**

**6m2 booths** are numbered:

- 1 – 4
- 10 – 11
- 14 – 21
- 28 – 30a
- 39 – 43

See room measurements below



- Booth 1: R-Biopharm
- Booth 2: aprimeo
- Booth 3: bioMérieux
- Booth 4: Mediphos
- Booth 5: Biocartis
- Booth 6: Biogio
- Booth 7: Qiagen
- Booth 8: ArcticZymes Technologies
- Booth 9: Seegene
- Booth 10: Copan
- Booth 11: BD Diagnostics
- Booth 12: Roche
- Booth 13: Roche
- Booth 14: LaCAR MDX Technologies
- Booth 15: MolGen
- Booth 16: SpeedX
- Booth 17: Labhelp Robotics
- Booth 18: LumiraDx
- Booth 19: Clean NA
- Booth 20: Medical Wire
- Booth 21: Altona
- Booth 22: Hologic
- Booth 23: Menarini Diagnostics
- Booth 24: Abbott
- Booth 25: NimaGen
- Booth 26: Goffin Molecular Technologies
- Booth 27: SD Biosensor
- Booth 28: Thermo Fisher
- Booth 29: Thermo Fisher
- Booth 30: LGC Bioscience Technologies
- Booth 30a: Genetic Signatures
- Booth 31: Eurogentec
- Booth 32: Medix Biochemica Group
- Booth 33: HiMedia / EWC Diagnostics
- Booth 34: GeneProof
- Booth 35: Wondfo Biotech
- Booth 36: UgenTec *Note: side parts should not be blocked!*
- Booth 37: Luminex a Diasorin company
- Booth 38: Elitech
- Booth 39: Axon Lab
- Booth 40: Molzym
- Booth 41: Magtivio
- Booth 42: Segers Medical Systems
- Booth 43: Cepheid

**Keizerzaal measurements:**

Square Meters:	400
General room measurements:	Width 21,23 / Length 19,90
Measurements lower section:	Width 21,23 / Length 7,60
Ceiling height higher part:	3,85
Ceiling height lower part:	2,70
Window surface:	corner to corner 5,58
Projection screen in the room:	Height 2,97 / Width 5,88
Door Height:	2,22 / Width 2,02

Pillar Keizerzaal to cabinet entrance kitchen: 2,01

Pillar Keizerzaal to wall Bar Maritime: 2,36

**Picke 1 measurements:**

Square Meters Pické 1:	374
Pické 1 L x W x H:	17 x 23 x 5
Ceiling heights:	6,48
Up to Chandelier:	4,67
Ceiling height rails dividing walls:	5,46
Ceiling height outer ring (wood):	3,38
Ceiling height outer ring:	3,52
Door Pické 1: Height 2,22 / Width 1,95 (between rails 1,72	
Walls landside per surface:	Width 7,42
Walls curves land side:	Width 6,35
Walls middle Pické 1:	Width 6,90
Walls with tapestry:	Width 4,69
Hanging points Max weight:	800 kg per hanging point
Diameter eye hanging point:	0.4
Measurements trusses:	Pické 1: 7,80 x 5,40

Measurements opening translator cabins:

Height 1,40 / Width 1,56 beam in the centre 0,68